PORTUGUÊS

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ENGLISH

FUNCTION AND INDICATION

The Chimiolux has as function to evidence bacterial plates or to evidence infiltration in teeth that have suffered any type of breaking. The product is indicated for evidence procedures in odontological clinics.

⚠ WARNING

This product is of restricted use of the dentists and professionals of the dentist area.

COMPOSITION

Water and methylene blue.

COMPARISONS BETWEEN THE MODELS

- **Chimiolux 5**
It has 0.005% of methylene blue concentration.
- **Chimiolux 10**
It has 0.01% of methylene blue concentration.

ACTION MECHANISM

The Chimiolux, when in contact with the tooth, is quickly absorbed, detaching the existing plates or breakings.

PRESENTATION FORMS

The Chimiolux is commercialized in the following Sets:

10 x 1,5 mL Set (pipe):

- 10 applicators (pipe type - as illustrated above) with 1,5 mL of Chimiolux 5 or Chimiolux 10, each;
- 1 User Manual.

10 x 1,0 mL Set (syringe):

- 10 applicators (syringe type - as illustrated above) with 1,0 mL of Chimiolux 5 or Chimiolux 10, each;
- 10 applicators peaks (as illustrated above);
- 1 User Manual.

STORAGE AND TRANSPORT

Keep in a place protected from direct light, heat and far from the children reach. Store the product in ambient temperature.

USE INSTRUCTIONS

- 1) Open the external packing and verify if the pipe packings or syringes are unbroken;
- 2) Open the pipe packing or syringe, only at the product application moment;
- 3) With the pipe a common syringe can be used to remove the product or to insert it in an applicator of local anesthetics (common in odontological offices). With the syringe, it is enough to place the applicator peak, which follows with the product;
- 4) Apply the product directly on the place where you want to identify the plate or infiltrations;
- 5) After 60 seconds, wash the application place with water and observe if the blue coloration is detached in the intended areas;
- 6) After the plate removal or infiltrations correction, the product must not remain next to teeth.

GENERAL INSTRUCTIONS OF SECURITY

- Restricted use of the Dentist Surgeon professional.
- The application of this product must be restricted to the treatment area;

- This product must be manipulated by the operator with protection gloves.

COLLATERAL EFFECT

There are no stories of collateral effect of this product.

ADVERSELY INDICATIONS

- Do not mix this substance with strong oxides as permanganate and oxygen water;
- Do not inject, nor use the Chimiolux in sanguineous vases.

CARES

Use the Chimiolux due the validity.

LIMITATIONS

There are no stories of limitations.

⚠ IMPORTANT INFORMATION

Solutions of methylene blue in the concentrations between 0,005 - 0.01%, when submitted to red light, with wave length between 620 - 680 nm, generate a singlet oxygen. This oxygen destructs the cellular membrane, generating a bactericidal action in consequence.

VALIDITY

The product has validity of two years.

USED SYMBOLS

-  Manufacture Date
-  Identify the manufacture lot
-  Product Validity Date
-  Manufactured by
-  Keep dry

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